

Product Name :
Cup Mill

Product Code :
TEST5146RFQLAB-0002



Description :

Cup Mill for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Cup Mill is designed for preparing powdered substance for X-ray fluorescence(XRF) Analysis. This vibration mill guarantees a high degree of reproducibility in the preparation of samples.

Technical Parameters

Rotation rate (rpm) 910

Vessel volume (cm³) 100

Feed grain (mm) 75

Dimension (m) 0.55x0.75x0.94

Net weight (kg) 230

Power supply AC380, 3?, 50Hz

Nominal power (kw) 1.1

3-Phase AC 380V

Permissible deviation $\pm 5\%$

Frequency 50Hz

Power 1.1kW

We are a leading manufacturers of Cup Mill testing equipment for a range of industries, including materials science testing lab. We offer full range of Cup Mill testing equipment, including Cup Mill testing kit, Cup Mill testing machine, Cup Mill testing tools, Cup Mill test equipment, Cup Mill test machines, Cup Mill for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission

is to provide reliable, accurate, and easy-to-use Cup Mill testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India