

Product Name :
Grinding Cup Mill

Product Code :
TEST5146RFQLAB-0003



Description :

Grinding Cup Mill

Technical Specification :

This grinding cup mill is used in laboratories to crush and grind sample particles sizes smaller than 6mm. After the crushing and grinding process, the output powder fineness is 0.1mm (120-200 mesh). The grinding mill features a streamlined mechanical structure with a tight cover and continuously accepts feeding samples.

It is among the most advanced crush and grinding equipment available on the market today.

Technical Parameters

Grinding Cup diameter ?175mm

Feeding particle size ?6mm

Discharge fineness 0.1mm

Capacity 25kg/hr

Motor power 1.5kw

Dimensions 750x550x400mm

We are a leading manufacturers of Grinding Cup Mill testing equipment for a range of industries, including materials science testing lab. We offer full range of Grinding Cup Mill testing equipment, including Grinding Cup Mill testing kit, Grinding Cup Mill testing machine, Grinding Cup Mill testing tools, Grinding Cup Mill test equipment, Grinding Cup Mill test machines, Grinding Cup Mill for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Grinding Cup Mill testing equipment that helps customers build better, safer, and more

sustainable infrastructure.



Testing lab instruments India