

Product Name :
Cement Mortar Mixer 5L

Product Code :
TEST5147RFQLAB-0001



Description :

Cement Mortar Mixer 5L for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Mortar mixer is designed to mix cement, water and sand to prepare cement paste or mortar for specimens that used in normal consistency tests and setting time tests.

Mixing blade width: 135mm

Mixing bowl capacity: 5L

Mixing blade speed: Rotation, Revolution

Low speed: 140 ± 5 r/min, 62 ± 5 r/min

High speed: 285 ± 10 r/min, 125 ± 10 r/min

Power supply: 220V AC, 50Hz, 370W

The shaft connecting screw between the stirring blade and the stirring leaf is M18x1.5.

The mixing pot volume is 5L and mixing pot wall thickness is 1.5mm.

The working gap between the stirring blade and the stirring pot is 3 ± 1 mm.

The motor is a vertical type split horsepower dual speed motor with a power of 0.55/0.37kw.

We are a leading manufacturers of Cement Mortar Mixer 5L testing equipment for a range of industries, including materials science testing lab. We offer full range of Cement Mortar Mixer 5L testing equipment, including Cement Mortar Mixer 5L testing kit, Cement Mortar Mixer 5L testing machine, Cement Mortar Mixer 5L testing tools, Cement Mortar Mixer 5L test equipment, Cement Mortar Mixer 5L test machines, Cement Mortar Mixer 5L for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction

industry. Our mission is to provide reliable, accurate, and easy-to-use Cement Mortar Mixer 5L testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India