

**Product Name :**  
Cement Paste Mixer

**Product Code :**  
TEST5147RFQLAB-0002



**Description :**

Cement Paste Mixer for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Cement mixer is used to mix cement and water to prepare homogeneous paste and specimen for normal consistency tests and setting time tests.

Speed: Revolution(r/min)

Low: 62±5

Stop

High: 125±10

Mixing bowl capacity: 2.5L

Power: 380V 50HZ 370W

Size: 472x280x466mm

Stirring blade width: 111±0.15mm

Shaft connecting screw for mixing leaf & blade: M16x1

Stirring pot: max. inner diameter 160±1mm, max. depth 139±2mm

Stirring pot wall thickness: 1mm

Working gap between stirring blade and stirring pot: 2±1mm

Outer size: 472x280x458mm

We are a leading manufacturers of Cement Paste Mixer testing equipment for a range of industries, including materials science testing lab. We offer full range of Cement Paste Mixer testing equipment, including Cement Paste Mixer testing kit, Cement Paste Mixer testing machine, Cement Paste Mixer testing tools, Cement Paste

---

Mixer test equipment, Cement Paste Mixer test machines, Cement Paste Mixer for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Cement Paste Mixer testing equipment that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India