



Email : sales@testinglabinstruments.com

Product Name :
Jolting Apparatus for Cement Mortar

Product Code :
TEST5147RFQLAB-0003



Description :

Jolting Apparatus for Cement Mortar

Technical Specification :

Jolting apparatus is used for the jolting-shaping of cement mortar trial-sample preparation.

Amplitude: 15mm

Vibrating frequency: 1Hz

Working cycle: 60

Distance between the center of table and shaft: 800mm

Power: 220V 50Hz 70W

We are a leading manufacturers of Jolting Apparatus for Cement Mortar testing equipment for a range of industries, including materials science testing lab. We offer full range of Jolting Apparatus for Cement Mortar testing equipment, including Jolting Apparatus for Cement Mortar testing kit, Jolting Apparatus for Cement Mortar testing machine, Jolting Apparatus for Cement Mortar testing tools, Jolting Apparatus for Cement Mortar test equipment, Jolting Apparatus for Cement Mortar test machines, Jolting Apparatus for Cement Mortar for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Jolting Apparatus for Cement Mortar testing equipment that helps customers build better, safer, and more sustainable infrastructure.

Testing lab instruments India

