

Product Name :
Cement Jaw Crushers

Product Code :
TEST5149RFQLAB-0001



Description :

Cement Jaw Crushers for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Jaw crushers are widely used to crush cement clinker, ore or rock specimens in laboratories.

Jaw plates are made of Manganese steels for wear resistance.

If feed particle size is larger than 60mm, please choose 100x150 jaw crushers.

If feed particle size is smaller than 60mm, please choose 125x100 jaw crushers.

Inlet size: 100x60mm

Feed particle size: ?45mm

Final particle size: 3~15mm

Capacity: 45~550kg/h

Motor power: 1.5kW

We are a leading manufacturers of Cement Jaw Crushers testing equipment for a range of industries, including materials science testing lab. We offer full range of Cement Jaw Crushers testing equipment, including Cement Jaw Crushers testing kit, Cement Jaw Crushers testing machine, Cement Jaw Crushers testing tools, Cement Jaw Crushers test equipment, Cement Jaw Crushers test machines, Cement Jaw Crushers for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Cement Jaw Crushers testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India