

Product Name :
Steel Bar Bending Cold Punch Set

Product Code :
TEST5153RFQLAB-0001



Description :

Steel Bar Bending Cold Punch Set for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Bending cold punch set is a special tool necessary for cold bending test of steel bar. It can be used together with universal testing machine or bending machine.

Technical Parameters

Number of punch: 32

Number of punch block: 4

Diameter of punch: 6~128mm

Scope of application: 6~32mm diameter

We are a leading manufacturers of Steel Bar Bending Cold Punch Set testing equipment for a range of industries, including materials science testing lab. We offer full range of Steel Bar Bending Cold Punch Set testing equipment, including Steel Bar Bending Cold Punch Set testing kit, Steel Bar Bending Cold Punch Set testing machine, Steel Bar Bending Cold Punch Set testing tools, Steel Bar Bending Cold Punch Set test equipment, Steel Bar Bending Cold Punch Set test machines, Steel Bar Bending Cold Punch Set for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Steel Bar Bending Cold Punch Set testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India