



Email : sales@testinglabinstruments.com

Product Name :
Grinding and Polishing Machine Touch Screen

Product Code :
TEST5158RFQLAB-0003



Description :

Grinding and Polishing Machine Touch Screen

Technical Specification :

Disc diameter: 300mm
Abrasive paper diameter: 300mm
Disc speed: Stepless speed 100~1000r/min
Disc rotational direction: Clockwise or counterclockwise
Disc electromotor: DC brushless motor, 220V, 1.2kW
Head electromotor: Stepper motor, 200W
Head speed: Stepless speed 20~120r/min
Timing set: 0~99min
Sample Holding Qty.: 6pcs
Sample holder: ?25mm, ?30mm, ?40mm
Specification: Custom specification is acceptable
Pressure method: Individual pressure or central pressure
Individual pressure force: 0~50N
Central pressure force: 0~160N
Control and display: 7 inch LCD touch screen
Power supply: 220V, 50Hz, 8A

We are a leading manufacturers of Grinding and Polishing Touch Screen testing equipment for a range of industries, including materials science testing lab. We offer full range of Grinding and Polishing Touch Screen

testing equipment, including Grinding and Polishing Touch Screen testing kit, Grinding and Polishing Touch Screen testing machine, Grinding and Polishing Touch Screen testing tools, Grinding and Polishing Touch Screen test equipment, Grinding and Polishing Touch Screen test machines, Grinding and Polishing Touch Screen for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Grinding and Polishing Touch Screen testing equipment that helps customers build better, safer, and more sustainable infrastructure.



Testing lab instruments India