

**Product Name :**  
Automatic Consolidometer

**Product Code :**  
TEST5165RFQLAB-0001



**Description :**

Automatic Consolidometer for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

This automatic consolidometer use special pneumatic pressure controller, with low noise, each unit sample consolidation pressure can be controlled and tested separately, testing result can be transferred via software system easily to computer.

Pressure Capacity: Medium pressure: 0kN - 4.8kN(0-1600KPa).

High pressure:0-9.6KN (0-3200KPa)

Load accuracy: Medium pressure:  $\pm 100\text{kPa}$ , relative error  $\pm 1.0\text{kPa}$

High pressure:  $\pm 100\text{kPa}$  - 1600kPa, relative error  $\pm 1.0\%$

Specimen size: 30cm250cm2

Pneumatic control range: 0MPa - 0.9MPa

Pressure transducer: 0MPa - 1.0MPa

Converter channel: 1 - 16

Power supply: 220VAC 50Hz

We are a leading manufacturers of Automatic Consolidometer testing equipment for a range of industries, including materials science testing lab. We offer full range of Automatic Consolidometer testing equipment, including Automatic Consolidometer testing kit, Automatic Consolidometer testing machine, Automatic Consolidometer testing tools, Automatic Consolidometer test equipment, Automatic Consolidometer test machines, Automatic Consolidometer for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use

---

Automatic Consolidometer testing equipment that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India