

**Product Name :**  
Consolidometer High Pressure

**Product Code :**  
TEST5165RFQLAB-0004



**Description :**

Consolidometer High Pressure

**Technical Specification :**

The Consolidation Apparatus is used for compression test of the soil, to determine the relationship between deformation and compression of the soil, to calculate the unit settlement ,compression index ,resilient index as well as consolidation coefficient of the soil etc.

Main Technical Specification:

Specimen size: 30cm<sup>2</sup>, 50cm<sup>2</sup>

Lever ration: 24:1 20:1

Loading method: Weight

Load pressure: 4000KPa(30cm<sup>2</sup>), 2000KPa(50cm<sup>2</sup>)

Pressure range: 12.5KPa ~ 4000KPa

Type: Duplex: Triplex

Cell: 2pcs: 3pcs

Outer size(L\*W\*H): 65\*45\*90cm: 65\*68\*90cm

Total weight(KG): 195(include weight): 290(include weight)

We are a leading manufacturers of Consolidometer High Pressure testing equipment for a range of industries, including materials science testing lab. We offer full range of Consolidometer High Pressure testing equipment, including Consolidometer High Pressure testing kit, Consolidometer High Pressure testing machine, Consolidometer High Pressure testing tools, Consolidometer High Pressure test equipment, Consolidometer High Pressure test machines, Consolidometer High Pressure for Lab Testing Instruments. We are suppliers of

---

materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Consolidometer High Pressure testing equipment that helps customers build better, safer, and more sustainable infrastructure.



## Testing lab instruments India