

**Product Name :**  
Consolidomete Light-Duty

**Product Code :**  
TEST5165RFQLAB-0005



**Description :**

Consolidomete Light-Duty

**Technical Specification :**

The test can be performed in order to determine the correlation/ratio between the applied pressure and the resulting deformation, or between vertical deformation and time during which a constant pressure is applied.

**Technical Parameters**

Load to be added: 12.5kPa-400kPa

Lever ratio : 1:20

Specimen area : 30cm<sup>2</sup> x 2cm(H), 50cm<sup>2</sup> x 2cm(H)

Instrument Dimensions: 500 x 170 x 380 mm (L x W x H)

We are a leading manufacturers of Consolidomete Light-Duty testing equipment for a range of industries, including materials science testing lab. We offer full range of Consolidomete Light-Duty testing equipment, including Consolidomete Light-Duty testing kit, Consolidomete Light-Duty testing machine, Consolidomete Light-Duty testing tools, Consolidomete Light-Duty test equipment, Consolidomete Light-Duty test machines, Consolidomete Light-Duty for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Consolidomete Light-Duty testing equipment that helps customers build better, safer, and more sustainable infrastructure.



# Testing lab instruments India