



Email : sales@testinglabinstruments.com

**Product Name :**  
Le Chatelier Mould

**Product Code :**  
TEST5056RFQLAB-0001



**Description :**

Le Chatelier Mould for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Used for determination of soundness of cement, and class-A and class-B limes. The apparatus consists of the following :-

One split cylinder 30 mm dia x 30 mm high made of spring brass with two indicator arms.

Two glass plates.

One lead weight approximately 100 gms.

We are a leading manufacturers of Le Chatelier Mould testing equipment for a range of industries, including materials science testing lab. We offer full range of Le Chatelier Mould testing equipment, including Le Chatelier Mould testing kit, Le Chatelier Mould testing machine, Le Chatelier Mould testing tools, Le Chatelier Mould test equipment, Le Chatelier Mould test machines, Le Chatelier Mould for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Le Chatelier Mould testing equipment that helps customers build better, safer, and more sustainable infrastructure.

---

# Testing lab instruments India

