

Product Name :
Automatic CBR/UCS/ITS/ Marshal Testing Machine

Product Code :
TEST5176RFQLAB-0003



Description :

Automatic CBR/UCS/ITS/ Marshal Testing Machine for Testing Lab Equipment Tender Suppliers India

Technical Specification :

AUTOMATIC CBR/UCS/ITS/ MARSHAL TESTING MACHINE

Features:

- Automatic Control
- Programming for speed and load rate
- Silent motor and gearbox design
- Not mechanical chains driven
- USB Connectivity for communication to PC, downloading of test log data and for uploading firmware updates
- Durable membrane touch button console

Technical Specifications:

- Mould /specimen diameter 152mm
- Platen travelling distance 80mm
- Capacity 200KN
- Speed rate (user defined) 1-55mm/min
- Load Rate (User defined) 40-150KN/min
- Pre-programmed tests CBR/ UCS /Marshall/ ITS
- Power Rating 230V 750 Watts

Accessories:

-
- 1) CBR Press Cover
 - 2) CBR Controller Software
 - 3) CBR Penetration weight
 - 4) CBR Penetration Piston
 - 5) UCS Breaking Head
 - 6) Marshall Breaking Head, ITT Break Head
 - 7) Laptop (CPU 1.8GHz Intel core i7-8565U 16GB RAM 13.3inch).



Testing lab instruments India