



Email : sales@testinglabinstruments.com

**Product Name :**  
Electric Laboratory Concrete Mixer

**Product Code :**  
TEST5176RFQLAB-0004



**Description :**

Electric Laboratory Concrete Mixer

**Technical Specification :**

Electric Laboratory Concrete Mixer

The machine is designed as a metallic body with mixing function that is performed by the paddles present inside the drum

The machine is suitable for mixing small quantities for lab purposes; it offers both wet and dry mixing of concrete and other materials.

Ease of accessing the mixture and also emptying after the task is done.

Our laboratory concrete pan mixer machine has been developed with industrial grade SS and fitted with highly durable and rust proof functional components that make it much durable and effective.

**Specifications:**

Type: Pan

Capacity: 300 liters

Mixing Capacity: 160-200 liters

Dimensions: 1250x1200x1300

Power: 5.5 kW

Discharge: Manually Controlled

Voltage: 220-230V 1 phase



# Testing lab instruments India