

Product Name :
High Capacity Sieve Shaker

Product Code :
TEST5176RFQLAB-0005



Description :

High Capacity Sieve Shaker

Technical Specification :

High Capacity Sieve Shaker

This high capacity sieve shaker is designed for sieving considerable quantities of any material.

The screen shaker accepts up to 30 litres of sample material.

Sturdily constructed, the machine can hold six screen tray sieves and dust pan. The screen tray sieves are 457mm wide x 660mm long x 75mm deep.

Supplied complete with the following screen trays

26.5mm

19.0mm 9.5mm

4.75mm

Dust Pan (receiver)

Capacity: 30litres (60-70Kg)

Power supply: 220-240V 50-60Hz

Tray dimensions: 457 x 660x75mm

Testing lab instruments India

