

**Product Name :**  
Standard Digital Viscometer

**Product Code :**  
TEST5176RFQLAB-0010



**Description :**

Standard Digital Viscometer

**Technical Specification :**

Standard Digital Viscometer

Used for determining the viscosity of cut back bitumen and road oil. The apparatus consists of a tank fitted with a thermostat, rheostat, an agitator, an immersion heater to take the water to the required temperature and a cooling coil for connection to the water supply.

Supplied with glass thermometer IP 8C range 0-45°C 0.2°C graduations

Graduated glass cylinder 100ml capacity

Go/no go gauge for 4mm diameter orifice

Go/no go gauge cup with diameter 10mm orifice

Cup with 10mm diameter orifice cup

Cup with 4mm diameter orifice

