



Email : sales@testinglabinstruments.com

**Product Name :**  
Balances

**Product Code :**  
TEST5176RFQLAB-0012



**Description :**

Balances

**Technical Specification :**

Balances

High Capacity, High Accuracy

Bright backlit LCD display - easy to read in all lighting conditions including bright sunlight.

Reads in multiple units g, lb, oz, ct, counting, % allows flexibility of weighing across many disciplines.

Filtering levels - can adjust balance response to differing environmental conditions for more controlled results.

Optional internal rechargeable battery - portability - can be used in various locations without access to a wall outlet.

Calibration is simple, ensuring that the balance will be calibrated.

Top Loading

a) Capacity 30000g

Resolution (g) 1

Pan size at least 600cm<sup>2</sup>

Mains and battery operated (batteries included)

Mains Voltage 220-240V 50-60Hz

---

# Testing lab instruments India

