

Product Name :
Balances 2000g

Product Code :
TEST5176RFQLAB-0013



Description :

Balances 2000g

Technical Specification :

Balances

Compact, lightweight body for easy portability.

Keypad calibration with external weight.

Back light for use in low-light conditions.

Rechargeable battery provides over eight hours of use on one charge.

Parts counting, percent weighing, HI/OK/LO comparator, remote zeroing, and full-scale tare.

Capacity: 2000g

Resolution (g): 1

Mains and battery operated (batteries included)

Mains Voltage: 220-240V 50-60Hz

