



Email : sales@testinglabinstruments.com

Product Name :
Water Bath For Marshal and Indirect Tensile Specimens

Product Code :
TEST5176RFQLAB-0016



Description :

Water Bath For Marshal and Indirect Tensile Specimens

Technical Specification :

Water Bath (For Marshal and Indirect Tensile Specimens)
Used to condition in water the Marshall ($60\pm 1^{\circ}\text{C}$) and other asphalt specimens.
Digital thermo regulator and temperature display, internal and external outer case in stainless steel. Complete with perforated base shelf and cover.

Capacity:

50-60Litres

20 -30litres

Digital thermo-regulator.

Internal and external outer case in stainless steel.

Continuous re-circulating system.

Complete with perforated base, shelf and cover.

Large, backlit LCD display to show set and actual temperature.

Safety:

Over-temperature/low liquid alarm

Redundant safety thermostat

Technical Specifications:

Recirculation YES

Temperature range: ambient to 600C
Accuracy $\pm 1^{\circ}\text{C}$
Resolution 0.1°C
Electrical 230v/50-60Hz).



Testing lab instruments India