

**Product Name :**  
Water Bath For Marshal and Indirect Tensile Specimens

**Product Code :**  
TEST5176RFQLAB-0016



**Description :**

Water Bath For Marshal and Indirect Tensile Specimens

**Technical Specification :**

Water Bath (For Marshal and Indirect Tensile Specimens)

Used to condition in water the Marshall ( $60 \pm 1^\circ\text{C}$ ) and other asphalt specimens.

Digital thermo regulator and temperature display, internal and external outer case in stainless steel. Complete with perforated base shelf and cover.

Capacity:

50-60Litres

20 -30litres

Digital thermo-regulator.

Internal and external outer case in stainless steel.

Continuous re-circulating system.

Complete with perforated base, shelf and cover.

Large, backlit LCD display to show set and actual temperature.

Safety:

Over-temperature/low liquid alarm

Redundant safety thermostat

Technical Specifications:

Recirculation YES

---

Temperature range: ambient to 600C  
Accuracy  $\pm 1^{\circ}\text{C}$   
Resolution  $0.1^{\circ}\text{C}$   
Electrical 230v/50-60Hz).



**Testing lab instruments India**