



Email : sales@testinglabinstruments.com

**Product Name :**  
ASTM Laboratory Thermometers

**Product Code :**  
TEST5176RFQLAB-0018



**Description :**

ASTM Laboratory Thermometers

**Technical Specification :**

| ASTM LABORATORY THERMOMETERS |          |               |           |
|------------------------------|----------|---------------|-----------|
| RANG                         | ASTM REF | RESOLUTION 0C | IMMERSION |
| +19 - +27oc                  | 17c      | 0.1           | total     |
| +18-+28oc                    | 24c      | 0.2           | 90mm      |
| -8 to +32oc                  | 63c      | 0.1           | total     |
| +34 to +42oC                 | 18c      | 0.1           | total     |
| +10 to +55oC                 | 76c      | 0.5           | 93mm      |
| -2 to 80oC                   | 15c      | 0.2           | total     |
| -5 to 110oC                  | 9c       | 0.5           | 57mm      |
| +90 to 170oC                 | 35c      | 0.2           | 50mm      |
| -1 to 175oC                  | 113c     | 0.5           | total     |
| Dial thermometers            |          |               |           |
| Range 0 to 260oC             |          |               |           |
| 50mm dial diameter           |          |               |           |

---

# Testing lab instruments India

