



Email : sales@testinglabinstruments.com

Product Name :
ASTM Laboratory Thermometers

Product Code :
TEST5176RFQLAB-0018



Description :

ASTM Laboratory Thermometers

Technical Specification :

ASTM LABORATORY THERMOMETERS			
RANG	ASTM REF	RESOLUTION 0C	IMMERSION
+19 - +27oc	17c	0.1	total
+18-+28oc	24c	0.2	90mm
-8 to +32oc	63c	0.1	total
+34 to +42oC	18c	0.1	total
+10 to +55oC	76c	0.5	93mm
-2 to 80oC	15c	0.2	total
-5 to 110oC	9c	0.5	57mm
+90 to 170oC	35c	0.2	50mm
-1 to 175oC	113c	0.5	total

Dial thermometers
Range 0 to 260oC
50mm dial diameter

Testing lab instruments India

