

**Product Name :**  
Automatic Water Distiller

**Product Code :**  
TEST5176RFQLAB-0019



**Description :**

Automatic Water Distiller

**Technical Specification :**

Automatic Water Distiller

All-in-one, automatically-controlled distillation water distiller made of stainless steel with an inbuilt storage tank.

**Specifications:**

Electronic control system (automatic water and power cut-off when tank is full and automatic start, Stainless steel heaters

Thermostatic low water cut off system

Original design, which allows to produce distilled water with temp. about 35-40° C without separate chiller

Electronic leakage protection system with solenoid valve

Easy-access evaporator, condenser-tank for effortless maintenance (cleaning, etc)

Easily replaceable condenser

Automatic level control in the evaporator

Visual control of the distillate level in the tank (through glass level-tube on the front panel)

All parts that contact water or steam are made of high-quality stainless steel

Stainless steel tubular heater with extended resource

Thermostatic low-water cut-off safety system.

Output: 4 liters/ hour

Water consumption better: 70 liters/ hour or

---

Storage tank Capacity tank: 8 liters  
Distillate conductivity: 2.5 $\mu$ S/m  
Power Supply: 230V 50-60Hz, single phase



**Testing lab instruments India**