

Product Name :
QUV Accelerated Weathering Tester

Product Code :
TEST5176RFQLAB-0022



Description :

QUV Accelerated Weathering Tester

Technical Specification :

QUV Accelerated Weathering Tester
Specimen Orientation 75o from the horizontal
UV florescent lamps (T12x121 cm long) @ 40watts each
Onboard radiance sensor Broadband(300-800nm)
Black panel temperature oC
Light cycle temperature
Condensation cycle temperature: 35 - 80
Specimen Exposure Area 20x50 (2 x each side) 4000cm2 total
Specimen capacity 48 specimens (75x150mm)
Inlet water pressure 0.2-5.5 bar
Inlet water purity Tap water
Water consumption 5 litres per day
Electrical requirement 230V±10%

Accessories:

Sample holders
Standard flat panel holder 75x150mm
Medium Flat Panel Holder 100x150mm
Large flat panel Holder 150mm x 150mm



Testing lab instruments India