

Product Name :
Automatic Digital Shear Box

Product Code :
TEST5176RFQLAB-0035



Description :

Automatic Digital Shear Box

Technical Specification :

Automatic Digital Shear Box
Suitable for samples containing particles up to 20mm largest dimension
To test sample sizes 150mm and 300mm square
Vertical pressure applied by an automatically programmable hydraulic system.
Shear vertical force 100KN
Speed range 0.00001 to 11.0000 mm/min
Up to 50 consolidation steps
Travel 75mm

Accessories

150mm square sample insert for 300mm shear box

Data processing software:

- i) Direct and residual shear Geo-Analysis template
- ii) Direct Shear Geo Analysis template.

Testing lab instruments India

