

**Product Name :**  
Automatic Digital Shear Box

**Product Code :**  
TEST5176RFQLAB-0035



**Description :**

Automatic Digital Shear Box

**Technical Specification :**

Automatic Digital Shear Box  
Suitable for samples containing particles up to 20mm largest dimension  
To test sample sizes 150mm and 300mm square  
Vertical pressure applied by an automatically programmable hydraulic system.  
Shear vertical force 100KN  
Speed range 0.00001 to 11.0000 mm/min  
Up to 50 consolidation steps  
Travel 75mm

**Accessories**

150mm square sample insert for 300mm shear box

Data processing software:

- i) Direct and residual shear Geo-Analysis template
- ii) Direct Shear Geo Analysis template.

---

# Testing lab instruments India

