

Product Name :
Steel Transverse Flexural Accessory For Steel Tester

Product Code :
TEST5176RFQLAB-0042



Description :

Steel Transverse Flexural Accessory For Steel Tester

Technical Specification :

Steel Transverse Flexural Accessory For Steel Tester
Type Of Testing Machines: Flexure Testing Machine
Material: Mild Steel
Finishing: Polished
Size: 4X6feet
Automation Grade: Automatic
Frequency: 50 Hz
Surface Finish: Polished
Maximum Load (KN): 200
Maximum distance between bearer (mm): 800
Bearer length (mm): 120
Bearer diameter (mm): 50.

