

**Product Name :**  
Laboratory Ovens

**Product Code :**  
TEST5176RFQLAB-0046



**Description :**

Laboratory Ovens

**Technical Specification :**

Temperature control and display with digital thermo-regulator/indicator.

Mechanical safety thermostat to avoid over temperature.

High temperature uniformity and precision

Over temperature protection

Stainless steel internal and external lining "linen patterned"

Forced convection airflow

Supplied with two extractable grid shelves

Voltage: 230 V, 50-60Hz 1 Phase

1) Nominal Capacity 100 liters

Temperature oC 200

Air Convection: Forced

Number of Shelves: 2

2) Nominal Capacity: 220 liters

Temperature oC 200

Air Convection Forced

Number of Shelves 2

3) Nominal Capacity 400 liters

Maximum temperature 0C 200

Air Convection Forced

---

Number of Shelves 3  
Overheat protection  
Digital temperature control system.



**Testing lab instruments India**