

Product Name :
Aggreagate Abrasion Value Dorry Abrasion Machine

Product Code :
TEST5177RFQLAB-0001



Description :

Aggreagate Abrasion Value Dorry Abrasion Machine for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Aggreagate Abrasion Value - Dorry Abrasion Machine

It consists of a disc rotating about a shaft connected to a reduction gear box coupled to a motor. The disc rotates at 28-30 RPM.

Under the rotating disc is a tray with an outlet to facilitate the removal of sand. Two Conical Hoppers are mounted on a bracket fixed to the circular tray.

An arrangement is made for start and stop the flow of sand. Two containers with weights are supplied to keep the specimens pressed against the rotating disc.

Specifications:

600mm diameter cast iron grinding disc

Grinding disc rotation speed 20-30 revolutions per minute

Specimen moulds 2

Two trays

Two flat plates

Weights and clamps

Voltage 230V 50/60Hz



Testing lab instruments India