



Email : sales@testinglabinstruments.com

**Product Name :**  
Benchtop Burst Strength Tester

**Product Code :**  
TEST5178RFQLAB-0005



**Description :**

Benchtop Burst Strength Tester

**Technical Specification :**

Benchtop Burst Strength Tester

Suitable for testing the seal integrity of liquid, gas filled - pouches / blister packs. Leak and Seal Strength Tester, applicable to quantitative determination of seal strength, seal quality, burst pressure, seal integrity, compression resistance, torsion force and joint/disengaging force of various flexible packages, aseptic packages, plastic tamper- evident closures, flexible tubes, caps and other materials.

Multiple Test Modes & Automatic Testing Process

Positive pressure method.

Multiple test modes including burst test, creep test, creep-to-failure test and pressure retaining test.

Dual test methods of restraint distension and unrestraint distension at customer's choice.

Real-time display of test curves and automatic statistics of test data.

Test range can be selected for non-standard tests.

Main components and parts are supplied by well-known global manufacturers, with guaranteed performance.

Innovative Real-time Monitoring of Gas Inlet Flow

Real-time monitoring of gas inlet flow under high pressure ensures the uniform testing condition during testing process.

High repeatability of test data.



# Testing lab instruments India