

Product Name :
Mechanical Rotating Sample Divider Motorised

Product Code :
TEST5067RFQLAB-0001



Description :

Mechanical Rotating Sample Divider Motorised for Testing Lab Equipment Tender Suppliers India

Technical Specification :

For obtaining representative sample for sieve analysis. It is suitable for dry sands, fine powders, ores, refractory materials, mineral aggregates, pigments, powdered coal, soap, cement etc. Up to a maximum particle size of 3 mm.

The sampler consists of the following :-

A hopper with an adjustable orifice having six holes of 1.6 mm, 3 mm, 5 mm, 6 mm, 8 mm and 10 mm dia and a shutter.

A rotating spreader turning at approximately 440 rpm to distribute the sample over the trough. Suitable for operation on 220v, 50 Hz, single phase supply