



Email : sales@testinglabinstruments.com

**Product Name :**  
Capillary Tube Cutter

**Product Code :**  
TEST5189RFQLAB-0003



**Description :**

Capillary Tube Cutter

**Technical Specification :**

Capillary tube cutter  
Cutting all size of capillary tube without collapsing tubes inside diameter.  
For cutting capillary tubes without collapsing or narrowing the tube. Cuts all sizes of capillary tube.  
Hardened steel construction.  
Cuts all sizes of capillary tubes.  
Won't collapse tube.  
Rubber grip handle.



**Testing lab instruments India**