

Product Name :
Sieve Shaker

Product Code :
TEST5194RFQLAB-0009



Description :

Sieve Shaker

Technical Specification :

Sieve Shaker

Sieve Shaker is ideally used as a particle size analyser for governing the separation of particles according to their sizes. The instrument can handle pharmaceutical powders and granules as well as chemical processing materials, dyes, chemicals and is also used by the food processing industry. It is also used to test sand, cement and gravel.

Sieving range from 0.005 mm to 4 mm

Pressure of 20 Kpa

Time display:1-99 minutes

Sieve capacity:one 200 mm diameter

Power rating:1200 W

Voltage:230 V,50-60 Hz,1ph

Total Wieght: 25 Kg

Accessories of Sieve shaker

a)Noise reduction cabinet

Made from sheet steel and lined

Internall with soundproofing material

Overall Dimensions:870*672*1562 mm

Weight: 90 Kg

b) Steel base plate for floor-mounting.



Testing lab instruments India