

**Product Name :**  
Casagrande Apparatus

**Product Code :**  
TEST5194RFQLAB-0012



**Description :**

Casagrande Apparatus

**Technical Specification :**

Casagrande apparatus

Used to determine the moisture content at which clay soils pass from a plastic to a liquid state.

Comprises: removable brass cup, adjustable crank and cam mechanism, blow counter, and base. Different versions are available conforming to the various Standards in use.

They are identical in shape and differ, generally, from the type of base and weight of the cup.

Motorized versions :4Kg

Accessories

Metal grooving tool

Spare brass cup