

Product Name :
Cone Penetrometer

Product Code :
TEST5194RFQLAB-0013



Description :

Cone Penetrometer

Technical Specification :

Cone Penetrometer

Semi-automatic digital liquid limit penetrometer with vertical micrometric adjustment and electronic release mechanism

Used to determine the moisture content at which clay soils pass from a plastic to a liquid state. Two Version are available. One the standard one the other is semi-automatic with timer where the cone is allowed to freefall for a period of 5 seconds including one each of 50 gm. 100 gm. Weight, one penetration cone, preset counter measuring cup.

230V,50-60Hz,1ph