

Product Name :
Compaction Test

Product Code :
TEST5194RFQLAB-0015



Description :

Compaction Test

Technical Specification :

Compaction Test
Standard proctor mould & rammer
Complete mould with base and collar
Mould volume 944 Cm³
Internal diameter 101.6 mm
Body height 116.4 mm & weight 7 Kg
Rammer with diameter 50.8 mm, weight 2.49 Kg
Accessories
Universal Extruder: made of steel, capacity 50 KN
Ram & screw travel 265 mm & weight 25 Kg

