

Product Name :
Permeability Test

Product Code :
TEST5194RFQLAB-0022



Description :

Permeability Test

Technical Specification :

Permeability test constant head permeability cell
114 mm internal diameter ,12 (blaked) take off points

Transparent cell body

Anodised aluminium upper & lower plate

Weight 6.5 Kg

Accessories

Manometer stand

Manometer stand with 3 manometer tube

The metal stand comprises 3 constant bore tubes, a meter scale, nipples & connecting tubing,

Dimension :1150*200*50 mm

Weight :3 Kg

Constant level tank

Made acrylic glass, complete with inlet, outlet, overflow, connecting tubing for the cell and Attachment for wall mounting.

Dimension :250*150*160 mm

Weight: 2.6 Kg

Tamping rod

8 mm diameter

300 mm height



Testing lab instruments India