

**Product Name :**  
Digital Hot Plate

**Product Code :**  
TEST5194RFQLAB-0034



**Description :**

Digital Hot Plate

**Technical Specification :**

Digital hotplate

Body is fabricated out of thick mild steel duly finished with powder coating and SS Top. Heater made of Kanthal wire to operate on 220 V AC 50 Hz single Phase. Heat is controlled by Microprocessor Digital Temperature indicator cum controller.

Power supply: 220-240VAC ,50Hz,1ph

Diam.: 520x130mm

Construction: Robust metal with anodised aluminium plate, heated surface

Max, temperature (C): 300x1

Rated power: 1500 W