

Product Name :
Sieve Analysis Shaker Digital With All Accessory

Product Code :
TEST5194RFQLAB-0035



Description :

Sieve Analysis Shaker Digital With All Accessory

Technical Specification :

Sieve analysis shaker digital with all accessory

Sieve Shaker is ideally used as a particle size analyser for governing the separation of particles according to their sizes. The instrument can handle pharmaceutical powders and granules as well as chemical processing materials, dyes, chemicals and is also used by the food processing industry. It is also used to test sand, cement and gravel.

Dimension: 380x 440x1085mm

Power supply: 220- 240 VAC, 60Hz,1ph