

Product Name :
Ostwalds Viscometer

Product Code :
TEST5195RFQLAB-0001



Description :

Ostwalds Viscometer for Testing Lab Equipment Tender Suppliers India

Technical Specification :

The relationship between the intrinsic viscosity of a polymer solution and the relative molecular mass of the polymer allows this apparatus to be used for quantitative determination of relative molecular mass.

Comprises glass U-tube with two bulbs, connected by a capillary tube.

Size: 120x1 mm (length x internal diameter)

Overall height: 237 mm

Material: Glass.