

**Product Name :**  
Ostwalds Viscometer

**Product Code :**  
TEST5195RFQLAB-0001



**Description :**

Ostwalds Viscometer for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

The relationship between the intrinsic viscosity of a polymer solution and the relative molecular mass of the polymer allows this apparatus to be used for quantitative determination of relative molecular mass.

Comprises glass U-tube with two bulbs, connected by a capillary tube.

Size: 120x1 mm (length x internal diameter)

Overall height: 237 mm

Material: Glass.