

**Product Name :**  
Heating Mantles

**Product Code :**  
TEST5195RFQLAB-0002



**Description :**

Heating Mantles for Testing Lab Equipment Tender Suppliers India

**Technical Specification :**

Built-in adjustable heater control up to an element temperature of 450°C  
Magnetic stirring control up to 520 rpm with autoreverse facility  
Automatic stirrer bar trap; the stirrer unit has an independent power supply.  
Polypropylene outer case is resilient and chemically-resistant.  
Even at full power output, the exterior remains "cool-to touch".  
All heating mantles have support clamp for rods of up to 12 mm diameter, and are double fused, with earth (ground).