

Product Name :
Heating Mantles

Product Code :
TEST5195RFQLAB-0002



Description :

Heating Mantles for Testing Lab Equipment Tender Suppliers India

Technical Specification :

Built-in adjustable heater control up to an element temperature of 450°C
Magnetic stirring control up to 520 rpm with autoreverse facility
Automatic stirrer bar trap; the stirrer unit has an independent power supply.
Polypropylene outer case is resilient and chemically-resistant.
Even at full power output, the exterior remains "cool-to touch".
All heating mantles have support clamp for rods of up to 12 mm diameter, and are double fused, with earth (ground).