



Email : sales@testinglabinstruments.com

Product Name :
Standard Sieve Set

Product Code :
TEST5195RFQLAB-0003



Description :

Standard Sieve Set

Technical Specification :

Sieves are of the highest quality to ensure consistent fit, accurate specifications, and durable brass construction.
Sieve Pans nest underneath sieves to collect fines
Sieve Covers prevent sample loss during agitation
Verification Services provide NIST traceable documentation of opening sizes
Sieve no: 10, 22, 36, 44, 65, 80, 100, 120



Testing lab instruments India