



Email : sales@testinglabinstruments.com

**Product Name :**  
Standard Sieve Set

**Product Code :**  
TEST5195RFQLAB-0003



**Description :**

Standard Sieve Set

**Technical Specification :**

Sieves are of the highest quality to ensure consistent fit, accurate specifications, and durable brass construction.  
Sieve Pans nest underneath sieves to collect fines  
Sieve Covers prevent sample loss during agitation  
Verification Services provide NIST traceable documentation of opening sizes  
Sieve no: 10, 22, 36, 44, 65, 80, 100, 120



**Testing lab instruments India**