



Email : sales@testinglabinstruments.com

Product Name :
Porcelain Evaporating Dish with Pour Spout and Watch Glass

Product Code :
TEST5195RFQLAB-0004



Description :

Porcelain Evaporating Dish with Pour Spout and Watch Glass

Technical Specification :

Evaporating dishes (or crystallizing dishes) are shallow bowls with a flat bottom and pour spout that may be used as reaction vessels, or for the separation of the solute from a solution through crystallization.

They are made from a range of thermally stable materials, but the ones at RHS are almost always porcelain.

A watch glass is a simple concavo-convex glass circle that is used to cover reaction vessels to prevent spattering during chemical reactions or to cover larger storage containers (e.g., beakers, evaporating dishes) temporarily to reduce loss of chemicals from evaporation.

Simple concavoconvex glass circle

Transparent

Capacity (mL): 3,250

Diameter (mm): 305

Height (mm): 95

Watch Glass, 150mm,

Fits 2000mL Beaker

Testing lab instruments India

