



Email : sales@testinglabinstruments.com

Product Name :
Weighing Balance

Product Code :
TEST5195RFQLAB-0005



Description :

Weighing Balance

Technical Specification :

Analytical electronic chemical balance for industrial and scientific field.
Sensor's Technology method: Strain gauge
Display: LCD
Wind proof (shield)
Dimension: 170*190*170mm
Outlook size: 270*195*285mm
Operation Key: ON/OFF, COUNT, ENTER, CAL, UNIT



Testing lab instruments India