


Product Name :
Pipette Filler Rubber Bulb (2 Sizes)

Product Code :
TEST5195RFQLAB-0014



Description :

Pipette Filler Rubber Bulb (2 Sizes)

Technical Specification :

It's efficient for one-handed operation.meet your needs perfectly well,such as control air evacuation, uptake and dispense liquids,etc

2pcs natural rubber made lab pipette filler bulbs with 3 Valves

Sizes of pipet bulb:5.4cm(2.1")in diameter, 13cm(5.1") in total length.

Accommodates pipettes up to 50ml.

The pipet filler is made of natural rubber, and has three valves that each containing a chemical-resistant stainless steel ball.