

Product Name :
Drop point of Grease Apparatus Gas Heated

Product Code :
TEST5079RFQLAB-0001



Description :

Drop point of Grease Apparatus Gas Heated for Testing Lab Equipment Tender Suppliers India

Technical Specification :

For determining the drop point of greases i.e. the temperature at which the grease under goes transition from the solid or semi-solid state to liquid state.

It consists of a brass cup 7.5 mm internal dia at wider end, a metal case and cylindrical metal sleeve, a glass boiling tube, beaker, stirrer and 7 mm dia steel ball.

We are a leading manufacturers of Drop point of Grease Apparatus Gas Heated testing equipment for a range of industries, including materials science testing lab. We offer full range of Drop point of Grease Apparatus Gas Heated testing equipment, including Drop point of Grease Apparatus Gas Heated testing kit, Drop point of Grease Apparatus Gas Heated testing machine, Drop point of Grease Apparatus Gas Heated testing tools, Drop point of Grease Apparatus Gas Heated test equipment, Drop point of Grease Apparatus Gas Heated test machines, Drop point of Grease Apparatus Gas Heated for Lab Testing Instruments. We are suppliers of materials testing tools and equipment for the construction industry. Our mission is to provide reliable, accurate, and easy-to-use Drop point of Grease Apparatus Gas Heated testing equipment that helps customers build better, safer, and more sustainable infrastructure.

Testing lab instruments India

