



Email : sales@testinglabinstruments.com

**Product Name :**  
Laboratory Balance

**Product Code :**  
TEST5195RFQLAB-0017



**Description :**

Laboratory Balance

**Technical Specification :**

Analytical balance features spacious weighing chamber, and the open-door clearance facilitates easy access to the weighing pan.

With this it is possible to use laboratory glassware of various sizes.

New aluminium base is a guarantee for weighing system stability.

Dual-CLICK system allows tool-free assembly and disassembly of the weighing chamber, and a special design of guides for weighing chamber panes (no rims) makes keeping the device clean convenient.

Innovative construction solutions increase balance durability and reliability in long-term operation.

Advanced analytical balance with internal adjustment system and memory

Readability [d]: 0.1mg

Advanced functionalities: dosing, density determination, peak hold, animal weighing, parts counting



**Testing lab instruments India**

